


Operating Manual of Fingertip Pulse oximeter

Model: JZK-303
JZK-303T



● Thank you for purchasing the Pulse Oximeter product produced by our company.
● In order for you to use this product safely and correctly, please read the operating manual before use.
● Please keep the operating manual properly.

VER: V1.0
Publication date: July 2020

Section 1

Safety

1.1 Instructions for the safe Operation and Use of the Fingertip Pulse Oximeter

Do not attempt to service the Fingertip Pulse Oximeter. Only qualified service personnel should attempt any needed internal servicing.

Do not use the pulse oximeter in situations where alarms are required. The device has no alarms.

Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status and correct alignment at least every 2 hours.

SpO₂ measurements may be adversely affected in the presence of high ambient light. Shield the sensor area (with a surgical towel, or direct sunlight, for example) if necessary.

The following reason will cause interference and defibrillators

- High-frequency electrosurgical.
- Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line.

The patient has hypotension, severe vasoconstriction, severe anemia or hypothermia.

The patient is in cardiac arrest or is in shock.

Fingernail polish or false fingernails may cause inaccurate SpO₂ readings.

1

1.2 Warnings

WARNING: EXPLOSION HAZARD

Do not use the Oximeter in a flammable atmosphere where concentrations of flammable anesthetics or other materials may occur.

Do not throw batteries in fire as this may cause them to explode.

Do not attempt to recharge normal dry-cell batteries. They may leak, and may cause a fire or even explode.

Do not use the Pulse Oximeter in an MRI or CT environment.

CAUTIONS

Caution

Keep the operating environment free of dust, vibrations, corrosive or flammable materials, and extremes of temperature and humidity.

Do not operate the unit if it is damp or wet because of condensation or spills. Avoid using the equipment immediately after moving it from a cold environment to a warm, humid location.

Never use sharp or pointed objects to operate the front-panel switches.

The batteries must be taken out from the battery compartment if the device will not be used for a long time.

The device shall only be used if the battery cover is closed.

The batteries must be properly disposed according to local regulation after their use.

2

1.3 Definitions and symbols

Symbol	Description
	Type BF Equipment
	Batch code
	Date of manufacture
	Information of manufacture, including name and address
	Temperature limitation
	When the end-user wishes to discard this product, it must be sent to separate collection facilities for recovery and recycling.
	The information you should know to protect patients and medical staff from possible injury.
	The information you should know to protect the equipment from possible damage.
	The important information you should know

3

Section 2

Introduction

2.1 General

This chapter provides a general description of the Fingertip Pulse Oximeter including:

- Brief device description
- Products features

2.2 Brief Device Description

The pulse oximeter, based on all digital technology, is intended for noninvasive spot-check measurement of functional oxygen saturation of arterial hemoglobin (SpO₂). Advanced DSP algorithm can reduce the influence of motion artifact and improve measurement accuracy of low perfusion.

The Oximeter can be used to measure human Hemoglobin Saturation and heart rate through finger. The product is suitable for use in family, hospital (including clinical use in internist/surgery, anesthesia, pediatrics, intensive care and etc.), Oxygen Bar, social medical organizations, physical care in sports and etc.

2.3 Product Features

- Lightweight for carrying and easy to use.
- Manually adjust the direction of interface.
- Color OLED/TFT display simultaneous display for testing value and plethysmogram.
- Low perfusion down to 0.3%.

4


Visual alarm function, real-time spotchecks.

- Low battery voltage indicator.
- Automatically switch off.
- Standard two AAA 1.5V Alkaline batteries support more than 20 hours continuous work.

Section 3

Installation, Setup and Operation

3.1 Description of the Front Panel




3.1.1 Parts of front and back panel
Table 3.1.1 Part Definition and Description

5

Item	Name	Description
1	Power button	Turn on the machine, direction change and parameter setting
2	OLED/TFT Panel	Display the SPO ₂ /PR data and Plethysmogram

3.2 Display

After switch on, the OLED/TFT display of the Oximeter is as follow



6

3.3 Parameter setting

Pressing the button (> 0.5s), oximetry enter the parameter setting, as shown below:

JZK-303	JZK-303T
Beep on	Alm on
SpO ₂ Lo 94	Pulse Beep on
PR Hi 130	SpO ₂ Lo 94
PR Lo 50	PR Hi 130
Restore off	PR Lo 50
+ / -	Restore off
Exit	Exit

Short press to enter the next setting, long press to adjust the parameters of the current cursor indicated item:
The parameters are described as follows:

Beep : Pulse oximeter alarm switch, ON is Enables / OFF is closed (JZK-303)
Alm : Pulse oximeter alarm switch, ON is Enables / OFF is closed (JZK-303T)
Pulse Beep : Pulse sound; ON is Enables/OFF is closed (JZK-303T)

SpO₂ Lo : Set lower limit of blood oxygen alarm

PR Hi : Set the upper limit of pulse alarm

PR Lo : Set lower limit of pulse alarm

Restore : Reset all settings to factory state

+ / - : Set the direction for digital adjustment. "*" Means the value increases upwards, "-" means the value decreases downwards

Exit : Exit the menu

7

3.4 Operation


3.4.1 Install batteries

Installing two AAA batteries into battery cassette in correct polarities and cover it.

WARNING: Do not attempt to recharge normal alkaline batteries, they may leak and may cause a fire or even an explode.

3.4.2 Turn the Pulse Oximeter on/off

Put one of fingers into rubber hole of the Oximeter (it is best to put the finger thoroughly with nail upward, then releasing the clamp).




Press power button to turn the Pulse Oximeter on. The oximeter will automatically be powered off when no finger in the device for longer than 8 seconds.

3.4.3 Read correspondent data from display screen

3.4.4 Display Description of OLED/TFT
The display interface of "OLED/TFT" can rotate four directions with four different display modes after pressing the power button for less than 0.5s. It is shown as below:

8



Note: When battery power is at lowest level, the battery capacity indicates symbol of "100%" in OLED/TFT, remind users of replacement of battery.

9

Section 4

4.1 Cleaning

Switch off the power and take out the batteries before cleaning. Keep the exterior surface (OLED/TFT display screen included) of the unit with a dry and soft cloth. Use 75% den sity of medical alcohol to clean the surface and use dry fabric with a little alcohol to avoid alcohol permeates into the device.

CAUTION: Don't use strong solvent. For example, acetone.
CAUTION: Never use an abrasive such as steel wool or metal polish.
CAUTION: Do not allow any liquid into the product, and do not immerse any parts of the device into any liquids.
CAUTION: Avoid pouring liquids on the device while cleaning.
CAUTION: Don't remain any cleaning solution on the surface of the device.

Section 5

Troubleshooting and Maintenance

5.1 Maintenance

- Replace the batteries timely when battery indication is low. Clean surface of the pulse Oximeter before it is used in diagnosis for patients.

10

- Remove the batteries inside the battery cassette if the Oximeter will not be operated for a long time.
- It is better to preserve the product in a place where ambient temperature is 10-40°C and humidity is 10%-80%.
- Regular inspection to make sure that no obvious damage existed to affect the safety and performance of the device.
- No flammable substance, overtop or lower temperature and humidity existed in operation conditions.

5.2 Troubleshooting

Problems	Possible Reason	Resolutions
Oxyhemoglobin or heart rate can not be shown normally	1. Finger is not plugged correctly. 2. Patient's perfusion is too low to be measured.	1. Retry by plugging the finger. 2. Try some more times. If you can make sure about no problem existing in the product. Please go to a hospital timely for exact diagnosis.
Oxyhemoglobin or heart rate is shown unstably	1. Finger might not be plugged deep enough. 2. Finger is trembling or patient's body is in movement status.	1. Retry by plugging the finger. 2. Try not to move, let the patient keep calm.

11

- 1. Power of batteries might be inadequate or not be there at all.
- 2. Batteries might be installed incorrectly.
- 3. The Oximeter might be damaged.

The Oximeter can not be powered on

- Please replace batteries.
- Please reinstall the batteries.
- Please contact with local customer service center.

The screen are suddenly off

- The product is automatically powered off when no signal is detected longer than 8 seconds.
- Power quantity of the batteries is exhausted.

- Normal.
- Replace the batteries.

Fingertip Pulse Oximeter

Specifications:

Physical Characteristics

Machine:

Dimension: 60 mm (L) x 35 mm (W) x 32 mm (H)

Weight: 50g (including 2 AAA alkaline batteries)

Outer box: 80 mm (L) x 60 mm (W) x 52 mm (H)

Gross weight: 70g

Outer carton: 510 mm (L) x 210 mm (W) x 235 mm (H)

Gross carton: 15kg

12

Classification

Anti-electronic Shock Type: Internally powered equipment

Anti-electronic Shock Degree: Type BF equipment

EMC: Type B class I

Mode of operation: Continuous Operation Enclosure Degree of ingress Protection: IPX4

Power

Internal	2 AAA 1.5V Alkaline batteries
Power consumption	Smaller than 30mA (Normal)

Temperature and humidity

Operating Temperature	5°C-40°C
Storage Temperature	-10°C-50°C
Relative Humidity	15%-80% non-condensing

Alarm default value

Parameter	Value
Blood oxygen saturation	Bottom limit: 94
Pulse rate	Upper limit: 130 Bottom limit: 50

13

Specifications

Parameters	Value
Measurement range of Blood oxygen saturation	45-100%
Measurement range of pulse	25-250 BPM
Resolution	Blood oxygen saturation 1% Pulse rate 1 BPM
Measure Accuracy	Blood oxygen saturation 2% (70%-100%) Unspecified (45%-69%) Pulse rate 2 BPM

14

Warranty Card

- This product is within one year's warranty under normal use of non-human factors of quality problems. (Gifts do not belong to the warranty.)
- The warranty date is subject to the invoice of the product, please show the invoices when maintained.
- During the warranty period, the repair is free under the following circumstances:
 - Do not follow the instructions for use caused by the fault.
 - Due to self-repair or modifying that resulting in failure.
 - Due to natural disasters, dropping, hitting or improper voltage.
 - Appearance is naturally polluted due to use.

Certificate

Product Name : Pulse Oximeter

Type : JZK-303 / JZK-303T

Inspection date : See body label

Examiner : _____

The product has been tested and proved to be qualified for delivery

